



Pediatric blanket series products



HB05-537
Small lower body shape



HB05-530
Small extended type

★ There are also pediatric blankets and postoperative blanket products.


Product features

- The unique arrangement of liquid drain ports ensures safe and efficient heating;
- The inner film can be used to cover the patient's body surface, which is helpful to improve the heating effect;
- Pediatric inflatable blanket is specially designed for young patients, without the need to choose special instruments and equipment;
- Lower body blanket and small extension blanket, sizes suitable for small patients of all ages from newborn to adolescent.


Specialty and cardiac surgical blanket series products

Product features


- The catheter design can guide the balanced distribution of heat in the core and peripheral parts of the human body;
- Effective body surface rewarming after cardiac surgery can obviously reduce the application of vasodilators and reduce the coagulation function after cardiac surgery;
- Designed by ethylene oxide sterilization, it is more suitable for advanced sterile surgery departments.




HB06-560
Catheter chamber type



HB06-570
Surgical path type



HB06-630
Catheter type



HB06-645
Cardiac surgery type





Shenzhen Medlink Medical Electronics Co., Ltd
Address: Floor 2, Floor 4, Floor 5, No. 2, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, Shenzhen
Tel: +86 18929727044
Website: www.med-linket-corp.com



* Disclaimer: All registered trademarks, product names, models, etc. displayed in the above content are owned by the original owners or original manufacturers. This article is only used to illustrate the compatibility of Midea's products and has no other intention! All the above information is for reference only, and should not be used as a work guide for medical institutions or related organizations. Otherwise, any consequences will have nothing to do with our company.

SP-106EN
V1.0 20241009

Disposable Sterile Warming Blanket




The disposable sterile warming blanket provided by Medlink belongs to the power inflatable, which meets the requirements of sensory control ability of the anesthesia operating room in the hospital. It can effectively solve the phenomenon of hypothermia of surgical patients, reduce the probability of chills during the recovery period, and shorten the recovery time of patients.


Medlink can provide 24 types of warming blankets, which can provide different options for clinical needs (such as preoperative, intraoperative, postoperative, cushion blankets). It is equipped with specialist warming blankets (such as cardiology, interventional catheters, pediatrics, lithotomy position, etc.) according to special needs to meet the warming needs of all patients.



Preoperative blanket series products



HB01-110
Outpatient patient with
soft sock type




HB01-111
Preoperative blanket
outpatient type


Product features

- Soft socks wearing design to better retain temperature and prevent body heat loss;
- It is conducive to the rewarming of patients after general anesthesia and provides favorable conditions for the smooth operation;
- Use a warming blanket before surgery to keep the patient in a warm environment and eliminate fear and tension.


Intraoperative blanket series products




HB03-522
Upper body shape



HB03-523
Upper body enlargement




HB03-525
Lower body shape




HB03-540
Trunk type

Product features

- The design of intraoperative blanket can provide uniform heat distribution;
- Flexible and elastic materials can meet the needs of various surgical positions;
- Non-inflatable foot pad protects heat-sensitive feet and calves from being burned;
- The attached transparent head towel can keep warm airflow around the intubated patient's head, and enable doctors to clearly observe the patient;
- Light weight and easy to handle after use.



HB03-542
Torso dual air inlet type




HB03-610
Systemic surgical type


Postoperative blanket series products

Product features


- After operation, the contact area of the blanket is larger, the edge is fully inflated, and the patient's body is fully insulated;
- Effectively shorten the time required for patients to wake up, reduce the incision infection rate and reduce postoperative complications;
- Inflation and heating have the highest efficiency, which can restore the patient's body temperature to normal in the shortest time.




HB02-300
Full body type



HB02-305
Thoracic path type




HB02-310
Pediatric type




HB02-315
Multipath type


Cushion blanket series products




HB05-555
Pediatric blanket
type




HB05-550
Pediatric enlarged
cushion blanket type




HB04-575
Spinal cushion
blanket type



HB04-585
Cushion blanket type
of lithotomy position



HB04-545
Adult cushion
blanket type



HB04-635
Full path mat
blanket type

Product features

- Before operation, lay the blanket on the operating table. It is convenient for rapid heating and saves preparation time;
- It is suitable for almost all types of surgery, and the unique design of the blanket will not cause any obstacle to the operation of medical staff;
- The novel design of the drainage hole can prevent the liquid from accumulating at the local compression point when the patient is lying on the blanket, and prevent the heating of the blood parts that may be deficient;
- The material is soft, X-ray transparent, and has no electromagnetic interference. The air outlet of the array ensures uniform heat transfer.