

Ordering Information

Order Code	Product Pictures	Descriptions	Packaging
Disposable EEG/EMG Electrodes			
B143-DP-10-60		DIN 1.5, Disposable Cup Electrode, 1.5 m	sets of 10, 1 x10 colors 10 colors: Black, White, Red, Yellow, Blue, Green, Purple, Orange, Brown, Beige
B143-DP-10-80		DIN 1.5, Disposable Cup Electrode, 2 m	
B143-DP-10-100		DIN 1.5, Disposable Cup Electrode, 2.5m	
SL0760A		Cup diameter: 9.2 mm, Cup height: 2.6 mm	20pcs/bag
Resuable EEG/EMG Electrodes			
B143-P-BLA-30		DIN 1.5, Cup Electrode, 0.75 m, Black	5pcs/bag
B143-P-BLA-60		DIN 1.5, Cup Electrode, 1.5 m, Black	
B143-P-GRN -30		DIN 1.5, Cup Electrode, 0.75 m, Green	
B143-P-GRN -60		DIN 1.5, Cup Electrode, 1.5 m, Green	
B143-P-RED-30		DIN 1.5, Cup Electrode, 0.75 m, Red	
B143-P-RED-60		DIN 1.5, Cup Electrode, 1.5 m, Red	
B143-P-10-48		DIN 1.5, Cup Electrode, 1.2m	sets of 10, 1 x10 colors 10 colors: Black, White, Red, Yellow, Blue, Green, Purple, Orange, Brown, Beige
B143-P-10-72		DIN 1.5, Cup Electrode, 1.8m	
B143-P-10-96		DIN 1.5, Cup Electrode, 2.4m	
B143-P-10-120		DIN 1.5, Cup Electrode, 3m	
B143-S-10-72		DIN 1.5, Snap Electrode, 1.8m	sets of 10, 1 x10 colors 10 colors: Black, White, Red, Yellow, Blue, Green, Purple, Orange, Brown, Beige
B143-S-10-96		DIN 1.5, Snap Electrode, 2.4m	
B143-S-10-120		DIN 1.5, Snap Electrode, 3m	
B143-S-5I-72		DIN 1.5, Snap Electrode, IEC, 1.8m	sets of 5, 1 x5 colors IEC colors: Black, Yellow, White, Red, Green, White AHA colors: Black, Brown, White, Red, Green, White
B143-S-5A -72		DIN 1.5, Snap Electrode, AHA, 1.8m	
B143-T-BLA-60		DIN 1.5, Alligator Clip Electrode, 1.5m, Black	5pcs/bag
B143-T-RED-60		DIN 1.5, Alligator Clip Electrode, 1.5m, Red	
CT0086A		EEG Electrode Cap Kit, Cu/AgCl	10pcs/bag
29-Lead EEG/EMG Electrode Set			
94061381B		Disposable EEG Electrode Set, 8 Leads (4+4), 0.5 m	1 set/bag
B113-DP-21-20		Disposable Cup Electrode Set, 21 Leads (11 + 10), 0.5 m	1 set/bag
B0070A		DIN 1.5, 29 Leads Adapter Cable for 94061381B and B113-DP-21-20, 2 m	1 set/bag

Declaration: Shenzhen Med-Link Electronics Tech Co., Ltd. ("MedLinket") hereby declares that all registered trademarks, product names, model numbers, etc., mentioned in this document belong to respective original owners or manufacturers. The information referenced herein is intended to clearly describe the compatibility features of MedLinket products for reference only and should not be construed as an operational guideline for medical institutions or related entities. If any adverse consequences cause from misunderstandings or incorrect use of the information, MedLinket shall not be held liable.



EEG/EMG Electrodes

Product Features

- Compatible with EEG, EMG, PSG, and neurophysiology systems for reliable bioelectrical signal transmission.
- Color-coded leads in multiple lengths and connector types allow easy identification, enhance clinical efficiency, and meet the needs of different application scenarios.
- Disposable and reusable cup electrodes are available to meet infection control standards and cost-effectiveness for long-term monitoring.
- High-quality gold plating and Ag/AgCl coating reduce electrode-skin impedance, promote stable signal acquisition, minimize artifacts, and enhance EEG signal quality.
- Constructed from medical-grade TPU with superior durability and flexibility for reliable performance in long-term use.

Application Scenarios

- Neurology
- Neurosurgery
- Neurocritical Care
- Sleep Medicine
- Epilepsy Monitoring Unit
- Clinical Research



**OEM & ODM Services Available,
Tailored to Your Needs**