

Premium Selection



Compatible Mindray non-invasive blood pressure hose

Product Advantage

- ★ Premium TPU material, bend-proof steady air sealing, durable;
- ★ Grey hose for adult & pediatric use, fast sorting by color marking;
- ★ Compliant biocompatibility, latex-free, cost-effective.

Scope of Application

The cuff end connector is used in combination with the blood pressure cuff. The equipment end connector is connected with the monitor to measure the patient's blood pressure.

Product Parameter

Compatible Models	Mindray BeneVision N series,BeneView T series,ePM series,iPM series iMEC series,uMEC series,MEC series etc.		
Brand	Medlinket	Order Code	YA106A06-10
Description	Single Tube Adult/Pediatric,10 ft (3m)	OEM#	6200-30-09688
Weight	110g/pcs	Shelf Life	5 years
Package	1pcs/bag	Warranty Period	6 months

*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: zhaosr@medlinket.com | store: <https://med-linket-corp.com> WhatsApp: +852 6467 3105